

Guidelines for handling and storage of Chemotechnique Diagnostics haptens

Storage:

The haptens should be stored in a refrigerator, (2-8) °C or a freezer, -18 °C (exceptional condition), as stated on the labels (long term storage). After use store the haptens according to recommended storage conditions.

Handling:

During transport, the haptens are subject to storage in room temperature, ≤25 °C, for normally up to 10 days and this does not harm the product. Patch testing conditions involves exposure to room temperature during the application of the haptens and this is considered normal handling.

Expiry:

Expiry dates are printed on the labels and we recommend renewing the product when the expiry period is ended.

Chemotechnique MB Diagnostics AB

**Quality Department
January 22th, 2018**